

A woman in white scrubs and glasses is shown in profile, using a handheld sprayer to clean a metal cart. The sprayer is connected to a black hose. The cart has a white tray on top and is on wheels. The background is a plain, light-colored wall.

**BECKER**

Product range:  
Endobox

**INFORMATION**

**3-2**  
**EN**

2018



Magnetic endoscope holder (EBS 8-C.51)



All the components of our endoscopy product range fit together perfectly. This makes work much easier: our practical Endoboxes and manoeuvrable, individually configurable equipment trolleys create mobility and ensure that instruments are transported safely and are properly placed ready for use. Magnetic holders make sure that endoscopes are stored safely and placed in the correct positions as needed.



Seals are always at hand to mark used and unused endoscopes.



*Can be machine-cleaned and disinfected, withstands temperatures up to 115°C*



**Endobox, grey**

500 x 535 mm, 133 mm high  
 W3086 70010

**Lid for Endobox, grey**

500 x 535 mm  
 W3086 70020

**Seal**, package containing 500 seals

red  W9060 71050  
 green  W9060 72050



131152 with Endoboxes/  
lids/seals



131300 with Endoboxes/  
lids/seals



**Endobox equipment trolley**

– suitable for pushing under other furniture –  
4 pairs of holding brackets for Endoboxes/Endobox-lids

4 swivel castors 75 mm ø,  
all with brakes;  
Dim.: 630 x 630 mm,  
height 790 mm  
□ 131150

4 swivel castors 125 mm ø,  
all with brakes;  
Dim.: 630 x 630 mm,  
height 855 mm  
□ 131152

**Endobox equipment trolley**

5 pairs of holding brackets for Endoboxes/Endobox-lids

4 swivel castors 75 mm ø,  
all with brakes;  
Dim.: 630 x 630 mm,  
height 955 mm  
□ 131160

4 swivel castors 125 mm ø,  
all with brakes;  
Dim.: 630 x 630 mm,  
height 1020 mm  
□ 131162

**Endobox transportation trolley**

8 pairs of holding brackets for Endoboxes/Endobox-lids

4 swivel castors 75 mm ø,  
all with brakes;  
Dim.: 630 x 630 mm,  
height 1560 mm  
□ 131180

4 swivel castors 125 mm ø,  
all with brakes;  
Dim.: 630 x 630 mm,  
height 1625 mm  
□ 131182


**Endobox transportation trolley, enclosed**

pair of hinged doors, lockable, 2 push handles,  
6 pairs of holding brackets for Endoboxes/Endobox-lids  
4 swivel castors 125 mm ø, 2 with brakes; wall protection strips  
[ 685 x 715 mm, height 1465 mm ] □ 131300

**Endobox transportation trolley, enclosed  
suitable for standard trolley washing plants** (tedesco CWA/WWA)

pair of hinged doors, lockable, 2 push handles,  
6 pairs of holding brackets for Endoboxes/Endobox-lids  
4 swivel castors 160 mm ø with CrNi-Steel casings,  
all with brakes; wall protection strips  
[ 685 x 715 mm, height 1505 mm ] □ 131301

*Material:*  
These trolleys are  
made of easy-to-clean  
and highly-hygienic  
chromium-nickel steel,  
AISI 304 (1.4301)



Pass-through  
cabinet



Endoscope  
reprocessing  
station, height-  
adjustable, with  
dosing device  
and Magnes  
accessories



For decades Becker quality steel furniture and equipment have been a fixed item in medical technology in clinics and hospitals.

Design, individual choice of colours as well as an extensive range of items and sizes allow you to make optimum use of your available space. The resulting concept is customized to your needs and logical through and through.



Metallwarenfabrik  
Walter H. Becker GmbH  
Anzenkirchener Straße 4  
84371 Triftern, Germany

Phone +49 (0)8562 9605-0  
Fax +49 (0)8562 9605-64  
info@becker-triftern.de  
www.becker-triftern.de

### Becker equipment:

- for the ward
- for the operating theatre
- for the central sterilisation unit
- for the endoscopy
- for all function rooms